



# Small RNA Library Prep Kit for NGS



## Product Introduction

This product is compatible for small RNA extraction from body fluids with a starting amount of 500 pg - 20 ng. It can finally be converted into a library that suitable for sequencing on Illumina and MGI platforms.

## Library Construction Workflow








## Advantages

Strategy	Problem
Library preparation for trace-amount small RNA (500pg)	Low concentration
Patented adaptor/linker modification and dimer cleavage technology	Linker dimer formation
UMI sequence in adaptor	PCR preference
Gel free	Lengthy gel excision
UDI adaptor	Easy to be cross-contaminated during sequencing
Comptible for sequencing on Illumina and MGI platforms	Single sequencing platform

## Product Characteristics

Single and sharp peak





## Order Information

Product Cat. No.	Product Name	Product Specifications
KD22101	Small RNA Library Prep Kit for NGS(MGI Platforms)	32 rxns
KD24101	Small RNA Library Prep Kit for NGS(Illumina Platforms)	32 rxns

## Contact Information

**Tel:** 0571-86700733

**Web:** <http://www.neohalobio.com>

**Support:** [support@neohalobio.com](mailto:support@neohalobio.com)

